
Circulate or not circulate, that is the question !

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
Why the liquid circulate Vapor bubbles from the vapor reservoir push up the liquid as these move upward.

Why bubbles emerge Due to the temperature difference between the top and lower ends, vapor pressure difference is created according to the Clausius-Clapeyron equation;

T_2 : higher temperature、 T_1 : lower temperature

$$\ln \frac{P_2}{P_1} = -\frac{L}{R} \left(\frac{1}{T_2} - \frac{1}{T_1} \right)$$

P_2 : vap. pressure at T_2 , P_1 : vap. pressure at T_1 , L : latent heat of vaporization, R : gas constant.

The liquid sealed in this heat-pipe is 80vol% ethanol aqueous solution. The figure at the right hand side shows that for this liquid the pressure difference of about 1/4 atm. exists when the temperatures of bottom and top parts are 20°C and 50°C respectively.

(cf : 760mmHg = 1 atmosphere.) Water boils at 100°C which means at this temperature, the equilibrium vapor pressure becomes 1 atmosphere.



Why bubbles emerge only from the vapor reservoir

When the lower part of the heat-pipe is immersed in hot water, the equilibrium vapor pressure of that part becomes high enough to generate a bubble (vapor phase) but because of the barrier of extra energy necessary to overcome the initial increase of free energy due to the creation of phase interface between liquid and vapor, vapor cannot nucleate without a vapor seed which only exists at the vapor reservoir.

$\Delta G = -\frac{4}{3}\pi r^3 G_v + 4\pi r^2 \sigma$

ΔG : total free energy change due to vaporization (must be negative in order for a bubble to grow). r : radius of a bubble.

G_v : free energy gain due to the phase change. σ : interface energy between liquid and vapor. r^* : the smallest diameter of the bubble that can expand (the critical radius for growth of a bubble).

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